

University of Texas Medical Branch Sleep Disorder Center Policy: 08.06.21 GRASS Keyboard F-Key Designations	Effective Date: January 2013 Revised Date: June 2012 Review Date: January 2013
--	--

## GRASS Keyboard F-Key Designations

- 
- Audience:** All personnel in the Sleep Disorder Center.
- 
- Purpose:** To assist the technologist with notations using the designated F-Keys on the keyboard.
- 
- Policy:** To use the F-Keys on the keyboard as a backup means for naming notations should the X-Keys key pad not function properly.
- 
- Procedure:** Keyboard F-Keys notations:

- F 2 – Lights Out
- F 4 – Eyes Open
- F 5 – Eyes Closed
- F 6 – Look Right/Look Left
- F 7 - Crying
- F 8 – Look Up/Look Down
- F 9 – Watching TV
- F 10 – Blink
- F 11 - Clench Jaw
- F 12 - Flex Right Foot
- SF 2 – Lights On
- SF 3 – Hold Breath
- SF 4 - Cough
- SF 5 – Movement
- SF 6 – Three Deep Breaths
- SF 7 – Snore
- SF 8 – Nasal Breath
- SF 9 – TIR/TOR
- SF 10 – Patient to Restroom
- SF 12 – Flex Left Foot

References:

GRASS Technologies, An Astro-Med Inc. Twin Recording and Analysis Software, Users Manual. Version 4.3. 2008

---

Shahzad Jokhio, M.D.  
Medical Director